## ChronoThane™ P 80A

Polyurethane (MDI)

CardioTech International, Inc.

## Message:

ChronoThane P is a family of aromatic ether based polyurethane elastomers. With a long history of reliable performance, this medical grade polymer has the versatility to be used across a broad range of applicational areas including catheters, ports and access devices.

These biocompatible materials possess characteristics such as low coefficient of friction, low extractables, dimensional stability, high impact resistance, and excellent tear strength.

ChronoThane P allows for ease of manufacturability and can be processed using conventional extrusion or injection molding equipment. These materials are available in hardnesses ranging from 75 Shore A to 75 Shore D.

AdvanSource Biomaterials synthesizes and manufactures medical grade materials offering the ability to tailor physical and mechanical characteristics to support and enhance your end product design.

These mechanical characteristic's, critical to the design and development of medical devices, can incorporate a wide range of physical and chemical properties while maintaining core characteristics such as biodurability and biocompatibility. In most materials, specialized characteristics such as the addition of colorant agents or antimicrobial properties (where applicable) can be added to the polymer to provide a homogenous material and limit secondary processing steps.

In addition, radiopaque agents may also be incorporated into the formula to provide additional product enhancements and may contain up to 40%, by weight, of a radiopaque agent thus allowing varied-scale visibility options.

With an expanding range of secondary operations including custom solution development, prototype coating capabilities, and project management services, ASB's expert team of chemists, scientists, engineers and industry professionals assist in every stage of customers' projects, from concept initiation through full-scale manufacture.

| General Information                   |                                 |          |             |  |  |
|---------------------------------------|---------------------------------|----------|-------------|--|--|
| Features                              | Aromatic                        |          |             |  |  |
|                                       | Biocompatible                   |          |             |  |  |
|                                       | Good Dimensional Stability      |          |             |  |  |
|                                       | Good Tear Strength              |          |             |  |  |
|                                       | High Impact Resistance          |          |             |  |  |
|                                       | Low Extractables                |          |             |  |  |
|                                       | Low Friction                    |          |             |  |  |
|                                       | No Animal Derived Components    |          |             |  |  |
|                                       |                                 |          |             |  |  |
| Uses                                  | Medical/Healthcare Applications |          |             |  |  |
| Agency Ratings                        | ISO 10993 Part 5                |          |             |  |  |
|                                       | USP Class VI                    |          |             |  |  |
| Forms                                 | Pellets                         |          |             |  |  |
| Physical                              | Nominal Value                   | Unit     | Test Method |  |  |
| Melt Mass-Flow Rate (MFR) (190°C/3.26 |                                 |          |             |  |  |
| kg)                                   | 2.0 to 26                       | g/10 min | ASTM D1238  |  |  |
| Water Absorption (Saturation)         | 1.2 to 1.4                      | %        | ASTM D570   |  |  |
| Hardness                              | Nominal Value                   | Unit     | Test Method |  |  |
| Durometer Hardness (Shore A)          | 80                              |          | ASTM D2240  |  |  |
| Mechanical                            | Nominal Value                   | Unit     | Test Method |  |  |
| Tensile Strength                      |                                 |          | ASTM D638   |  |  |
| Break                                 | 27.6 to 48.3                    | MPa      |             |  |  |

| 50% Strain                           | 4.48 to 5.86  | MPa  |           |
|--------------------------------------|---------------|------|-----------|
| 100% Strain                          | 5.52 to 6.89  | МРа  |           |
| 200% Strain                          | 8.27 to 9.65  | МРа  |           |
| 300% Strain                          | 11.4 to 13.8  | МРа  |           |
| Tensile Elongation (Break)           | 680 to 850    | %    | ASTM D638 |
| Injection                            | Nominal Value | Unit |           |
| Drying Temperature - Desiccant Dryer | 71.1 to 93.3  | °C   |           |
| Drying Time - Desiccant Dryer        | 3.0 to 4.0    | hr   |           |
| Drying time Desiceant Dryer          | 3.0 to 4.0    | 111  |           |
| Dew Point                            | -40.0         | °C   |           |

The information and data on this page are provided by manufacturers and document providers. SHANGHAI SUSHENG assumes no legal liability. It is strongly recommended to verify all technical data with material suppliers before final material selection. All rights belong to the original authors. If any infringement occurs, please contact us immediately.

## Recommended distributors for this material

## Susheng Import & Export Trading Co.,Ltd.

Tel: +86 21 5895 8519

Phone: +86 13424755533 Email: sales@su-jiao.com

No. 215, Lianhe North Road, Fengxian District, Shanghai, China

